

SKEYNDOR DERMALED® »» XPERT



ADVANCED
LED COSMETIC
PHOTOTHERAPY

SCIENCE AND LIGHT TO CARE FOR YOUR SKIN

SKEYNDOR DERMALED XPERT®

LED phototherapy cosmetic treatment with CE certification emitting 7 wavelengths to treat various skin conditions. A non-invasive technique that stimulates the natural processes of cellular regeneration.

Inspired by photosynthesis, this technology transforms light energy into visible skin benefits — pain-free, with no side effects or downtime, and without UV rays.

Its effectiveness is supported by clinical studies and approved by organisations such as the FDA for dermatological care.

DYNAMIC COSMETIC LED PHOTOTHERAPY

It uses 144 high-intensity, next-generation LEDs with specific wavelengths to stimulate regenerative metabolism at a cellular level.

An effective treatment that helps treat complex skin issues such as acne, post-acne scars, pigmentation, dullness and sensitive skin.



From the first treatment, the skin is revitalised, its texture and tone improved.

ADVANTAGES


Easily adapts to any skin type and causes no discomfort.

It is compatible with other cosmetic treatments and can be combined with radiofrequency, mesotherapy, oxygen therapy and post-procedure skin care (peelings, microneedling).

A versatile, high-quality device with visible results.


INTERACTION IN SKIN LAYERS

LED light phototherapy acts at different depths depending on its wavelength, modulating specific cellular processes in each layer of the skin.



SHORT WAVES 400–495 NM

Act on the **epidermis**, regulating bacterial and sebaceous activity with a purifying and soothing effect on the skin surface.



MEDIUM WAVES 495–600 NM

Penetrate the **superficial dermis**, modulating melanocyte activity, reducing hyperpigmentation and improving skin tone uniformity.




LONG WAVES 600–950+ NM

Reach the **deep dermis and hypodermis**, stimulating collagen and elastin synthesis, promoting tissue repair and skin firmness.

SKIN APPLICATIONS

The ideal treatment to address the **main skin conditions**.





SCIENCE AND LIGHT FOR OPTIMAL SKIN CARE

LED light acts at the cellular level, stimulating mitochondrial activity and promoting tissue regeneration.

Each colour penetrates different skin layers to target specific needs.

MAIN BENEFITS

- Stimulates collagen and elastin production.
- Reduces wrinkles and fine lines.
- Treats acne and pigmentation.
- Increases luminosity and vitality.
- Accelerates skin repair and healing.
- Soothes sensitive or irritated skin.
- Improves circulation and oxygenation.
- Repairs and protects damaged skin.
- No side effects or discomfort.

	COLOR	MAIN FUNCTION	SPECIFIC BENEFITS	SKIN TYPES/ CONCERNS
LONG WAVES	RED LIGHT (640NM) (epidermal action and reticular dermis)	Healing, regenerating, anti-inflammatory.	Stimulates collagen production. Improves skin firmness. Reduces wrinkles and fine lines. Decongests and promotes cellular repair.	Skin showing signs of ageing or in healing/recovery.
MEDIUM WAVES	YELLOW LIGHT (570NM) (superficial papillary dermis action)	Anti-inflammatory, anti-edematous, soothing, repairing.	Reduces redness and decongests the skin. Stimulates drainage and toxin removal. Improves capillary fragility.	Sensitive or reactive skin.
	GREEN LIGHT (523NM) (epidermal action)	Melano-regulating, antioxidant, repairing.	Enhances radiance and balances skin tone. Helps reduce skin hyperpigmentation.	Pigmented, uneven-toned skin.
SHORT WAVES	CYAN LIGHT (490NM) (epidermal action)	Superficial melano-regulating, antioxidant, sebum-regulating, bactericidal, hydrating.	It provides a dual benefit, balancing skin tone and radiance while also delivering a sebum-regulating and bactericidal action.	Sensitive skin with mild acneic lesions. Superficial post-inflammatory pigmentation.
	BLUE LIGHT (423NM) (epidermal action)	Sebum-regulating, bactericidal, antiseptic.	Reduces and balances excess sebum on the skin, acting on the bacteria responsible for acne breakouts (<i>C. acnes</i>).	Combination or oily skin with acne.
	VIOLET LIGHT (400NM) (epidermal action)	Sebum-regulating, bactericidal, oxygenating.	Combines the properties of red and blue light, stimulates the skin's regeneration processes, and reduces the proliferation of <i>C. acnes</i> bacteria.	Combination or oily skin with acne. Scars.
	WHITE LIGHT (400-700NM) (global action)	Global stimulation of cellular metabolism.	Stimulates global tissue regeneration, improves texture and boosts vitality.	All skin types for maintenance and revitalisation.

IMPROVE THE RESULTS OF YOUR FACIAL TREATMENTS*

BAGS AND WRINKLES
RED LIGHT




Volunteer 1 after 3 weeks of treatment.

ACNE
BLUE LIGHT



Volunteer2 after 2 weeks of treatment.

SPOTS
GREEN LIGHT



Volunteer 3 after 5 weeks of treatment.

ACNE & SCAR REDUCTION
PURPLE LIGHT



Volunteer after 6 weeks of treatment.

GENUINE REVIEWS*

The use of SKEYNDOR DERMALED XPERT® equipment has shown improvement in:

- Reduction in fine lines, puffiness, firmer, rejuvenated skin.
- Brighter skin, fewer wrinkles, reduced sunspots.
- Reduced acne scars and visible relaxing effect.
- Relaxing treatment with visible results.

*Photos and reviews provided by Beijing ADSS Development Co. Ltd. Study conditions not recorded.

MAIN REASONS TO CHOOSE SKEYNDOR DERMALED XPERT®

- **Latest-generation dynamic LED cosmetic technology**
Designed with high-intensity diodes for maximum clinical and aesthetic efficacy.
- **FDA-approved with international certification**
Safety and efficacy: meets the most demanding medical and aesthetic safety standards.
- **Therapeutic versatility**
Seven wavelengths for personalized treatment according to the skin's needs.
- **Non-invasive**
and 100% UV-free, painless, without side effects.
- **Clinically proven treatments**
Visible improvements in skin quality and care, with safe and measurable results.
- **Real results and proven effectiveness**
Visible improvements in acne, dark spots, wrinkles, and sagging from the first sessions.



- **Complements other professional treatments**
Enhances the facial results of peels, laser, radiofrequency, mesotherapy.
- **Easy and comfortable to use**
Automatic protocols and 20–30 min treatment sessions.
- **Personalized training**
Treatment protocols provided to optimize each treatment.



SKEYNDOR

SCIENCE creates BEAUTY®

SKEYNDOR HEADQUARTERS

Navas de Tolosa, 148
08223 Terrassa · Barcelona
Tel. +34 902 181 514

SKEYNDOR BARCELONA

Mallorca, 268
08037 Barcelona
Tel. +34 932 282 442

SKEYNDOR MADRID

Calle de Don Ramón de la Cruz, 17
28001 Madrid
Tel. +34 915 281 900

www.skeyndor.com



PYMEINNOVADORA



This material and all of the activities that SKEYNDOR carries out are fully aligned with the SDGs (Sustainable Development Goals).



Please be reminded that this catalogue is printed on recycled paper. If you no longer need it, please dispose of it in a blue recycling bin. Thank you very much.

8960CT/EN

